

## DATA SHEET

**Multiplex PCR Master Mix**

Master mix for multiplex PCR

Cat. Nº.	Amount
<input checked="" type="checkbox"/> POL-140XS	50 reactions
<input type="checkbox"/> POL-140S	100 reactions
<input type="checkbox"/> POL-140M	200 reactions
<input type="checkbox"/> POL-140L	500 reactions
<input type="checkbox"/> POL-140XL	1.000 reactions

**Shipping:**

Shipped on blue ice

**Storage Conditions:**

Store at -20 °C

**Additional Storage Conditions:**

Avoid freeze/thaw cycles.

**Shelf Life:**

12 months

**For *in vitro* use only!****Description:**

Multiplex PCR Master is specially designed for the set-up of multiplex PCR reactions. It contains an optimized composition of polymerase, nucleotides, MgCl<sub>2</sub> and stabilizing components in a specifically developed buffer system allowing the parallel amplification of a multitude of fragments in a single PCR assay. The master mix contains all reagents (except primer and template) in a 2x concentrated ready-to-use solution. The kit is recommended for use in clinical PCR reactions and highly suitable for multiple target gene amplification in a single tube. The high specificity and sensitivity of the mix is achieved by a chemically inhibited hot-start polymerase. Its activity is blocked at ambient temperature preventing the extension of nonspecifically annealed primers and primer-dimer formations at low temperatures during PCR setup.

**Kit contents:****2x Multiplex PCR Master (purple cap)**

Master mix of thermostable Hot Start DNA polymerase, dATP, dCTP, dGTP, dTTP, KCl, MgCl<sub>2</sub> and stabilizers.

**PCR grade water (white cap)****PCR Reaction Setup**

A reaction volume of 10-50 µl per assay is recommended for most PCR cyclers. Pipet with sterile filter tips and perform the setup in an area separate from DNA preparation or analysis. No-template controls should be included in all amplifications. Thaw, mix, and briefly centrifuge each component before use. Add the following components to a microcentrifuge tube:

**1. Prepare PCR master mix**

Note: Consider the volumes for all components listed next steps to determine the correct amount of water required to reach your final reaction volume.

Components	50 µL rxn	[ final ]
Water, grade PCR	To 50 µl	
2 X Multiplex PCR Master	25 µl	1X
Forward primer 1 (10 µM)	0,5 - 2,5 µl	0,1 – 0,5 µM
Reverse primer 1 (10 µM)	0,5 - 2,5 µl	0,1 – 0,5 µM
DNA template	x	10 pg – 1 µg**

\*\*genomic DNA: 1 ng-1µg; plasmidial ou viral DNA: 1 pg-1 ng

Cap each tube, mix, and briefly centrifuge the content.

## DATA SHEET

**Multiplex PCR Master Mix**

Master mix for multiplex PCR

**3. Incubate reactions in a thermal cycler.**

Recommended cycling conditions:

Step	Temp.	Time	
Initial denaturation	95 °C	10 min	
30-50 cycles	Denaturation	95 °C	15 - 30 sec
	Annealing <sup>1</sup>	58-64 °C	15 - 40 sec
	Elongation <sup>2</sup>	72 °C	1 min/kb
Final extension	72 °C	5 min	
Hold	4 - 8 °C		

1)The optimal annealing temperature (AT) can be calculated for each primer as following:  $AT = T_m - 5\text{ °C}$  with  $T_m = 2\text{ °C} \times (A+T) + 4\text{ °C} \times (G+C)$  Please note that primers should be designed to show minimal differences in there melting temperatures ( $T_m$ ).

2)The elongation time depends on the length of the fragments to be amplified. A time of 1 min/kb is recommended.

For optimal specificity and amplification an individual optimization of the recommended parameters may be necessary for each new template DNA and/or primer pair.