

**Calcium Chloride Solution 2.5 M**

Cat. No.	Amount
BU-103	100 ml

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at 4 °C

**Shelf Life:** 12 months

**Molecular Formula:** CaCl<sub>2</sub>

**Molecular Weight:** 110.98 g/mol

**CAS#:** 10043-52-4

**EC number:** 233-140-8

**Form:** liquid

**Concentration:** 2.5 M

**Description:**

Calcium Chloride, filtered sterile, is recommended for preparation of competent cells for transformation.

**Quality:**

- Molecular biology grade
- sterile filtered, filter pore size 0.22 µm
- DNase / RNase free

