

**Campylobacter jejuni Antigen**

native, Blood agar

Cat. No.	Amount
PR-BA139-S	250 µg
PR-BA139-L	1 mg

**For general laboratory use.****Shipping:** shipped on dry ice**Storage Conditions:** store at -80 °C**Additional Storage Conditions:** avoid freeze/thaw cycles**Shelf Life:** 12 months**Form:** liquid**Culture system:** Blood agar**Protein concentration:** Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.**Preservation:** None**Analysis:** SDS-PAGE and ELISA**Preparation:** Partially purified detergent extract of the membrane fraction. The antigen is presented in water. The main component is an outer membrane protein with approximately 45 kDa in SDS-PAGE.**Inactivation:** This product has been inactivated by disruption and detergent treatment. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.