

**Quantifoil™ R 1.2/1.3 on 200 gold mesh**

Cat. No.	Amount
X-101-Au200	100 pcs.

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a

**Applications:**

Support film for biological samples in Cryo-EM techniques

**Description:**

QUANTIFOIL™ R 1.2/1.3 carbon film has circular holes of 1.2 µm diameter with an interspace of 1.3 µm in both dimensions, applied on biocompatible 200 gold mesh. Suitable for magnifications around 50.000x and higher.